

## Modular Treatment Cabinet with One Bottom Drawer and Two Doors

### Model 7413M4

Treatment Cabinet has one bottom drawer and two doors with a shelf in the lower section, and two lockable doors with glass panes and adjustable glass shelf in the upper section. It comes equipped with a Kydex work surface and is the perfect complement to any exam room and treatment area.



**PRODUCTS ARE MADE IN THE USA**

#### Key features at a glance

- › 20 gauge steel construction
- › One adjustable glass shelf in upper section
- › Removable, non-adjustable shelf in bottom section
- › Lock on upper doors (with 2 keys)
- › Heavy duty concealed door hinges
- › Kydex top
- › Safety release latch on upper left door
- › Adjustable leveling glides
- › 3-year limited warranty

#### Specifications

Width: 25.25" (64.1cm)  
Height: 65" (165.1cm)  
Depth: 16.25" (41.3cm)  
Net Weight: 102lb (6.4kg)

Shipping Weight: 126lb (8.6kg)  
Shipping Dimensions:  
L28xW18xH71 (71.1cm x 45.7cm x 180.34cm)

#### Accessories

**7414M2** – Double door lock (factory installed)  
**7414M3** – Single drawer lock (factory installed)  
**7414M4** – Additional keys (x2)