

Manual Exam Table with Side Step and Four Storage Drawers

Model 7415M2

Manual Exam Tables with pneumatic-assist back have a 500lb weight capacity, and reversible side and front storage drawers to accommodate table placement in any room configuration. The table also features hidden stirrups, and an all-steel OSHA-compliant step with anti-slip guards. Ships assembled for easy installation.



Key features at a glance

- › 500lb capacity (226.8kg)
- › Removable, Seamless Upholstered Top with Prefixx® Protective Finish
- › Galvanized steel base, heavy gauge, powder coated
- › Comfort Grip Back Release Handle
- › One touch, easy-glide stirrups with 4 lateral positions and friction lock length adjustment
- › Slide out leg section with removable upholstered leg rest
- › Removable 9" by 9" plastic drip pan (22.9cm x 22.9cm)
- › Pull out steel non-slip side step 7" high (17.8cm)
- › 21" paper roll holder (53.3cm)
- › 2 front drawers on full extension roller bearing glides 18.25"x15.5"x4" each (46.4cm x 39.4cm x 10.2cm)
- › 2 reversible side storage drawers on full extension roller bearing glides 18.25"x15.5"x4" each (46.4cm x 39.4cm x 10.2cm)
- › 5-Year Limited Warranty

Specifications

- › Height: 32.75" (83.2cm)
- › Length: 58.5" – 74.5"
(148.6cm – 189.23cm)
- › Width: 28" (71.12cm)
- › Stirrups: Infinitely adjustable
- › Net Weight: 241lb (109kg)
- › Shipping Weight: 320lb (145kg)
- › Shipping Dimensions: L64xW30xH41
(163cm x 76cm x 104cm)

Accessories

- › **7416M1** IV Pole (Patient Left)
- › **7416M2** IV Pole (Patient Right)
- › **7416M3** Safety Rail Package
- › **7416M4** Removable Stainless Steel Drain Pan
- › **7416M5** Removable Stainless Steel Drain Pan w/Spout
- › **7416M6** 110v Drawer Warmer (Factory Installed)
- › **7416M7** 110v Outlet Package (Factory Installed)

Colors

							
Soft Linen	Smoky Cashmere	Warm Sand	Creamy Latte	Chocolate Truffle	Lemon Meringue	Cherry Blossom	Morning Fog
							
True Graphite	Classic Black	Fine Wine	Blue Skies	Lakeside Blue	Twilight Blue	Mint Leaf	Deep Forest



PRODUCTS ARE MADE IN THE USA

4767 Broadmoor Ave. SE, Suite #5, Kentwood, MI 49512

Phone: 616-514-1140 | Fax: 616-514-1141 | Toll Free: 877-514-1140 | www.medicus-health.com