



Fusion**ONE**TM

Series

Power Exam Tables

Accessible. Affordable.

ADA Access in Every Exam Room.






1 in 4

adults is living with
a disability.^[1]



1 in 2

over 65 have a
disability.^[2]



3 out of 4

are overweight
or obese.^[3]



The population over age
65 is projected to grow from
48 million today to

98 million
in 2060^[4].

Healthcare is changing

And traditional exam rooms are changing with it

Traditional box exam tables can be difficult to access, posing safety risks for patients and clinicians.

And budget constraints often make it difficult to place ADA-accessible power exam chairs in every exam room.

But not anymore!

With the ACCESSIBILITY of a power chair and the AFFORDABILITY of a box table, the FusionONE Power Exam table offers the features healthcare providers need most.


This hi-lo power exam table enables simple ADA access and ease of examination for patients of all ages and abilities.

The FusionONE Series of exam tables are available with a manual or power backrest and with the OneTouch Wheelbase™ System. With a variety of models and accessories, the FusionONE Series allows for easy healthcare facility standardization.



7418M2

[1] US CDC [2] Harris Family Center for Disability and Health Policy / Western University of the Health Sciences [3] Modern Healthcare [4] The future is not what it used to be: Thoughts on the shape of the next normal, Kevin Sneader and Shubham Singhal, McKinsey & Company April 2020DC.org



	7417M1	7417M2	7417M3	7417M4	7417M5	7418M1	7418M2	7418M3	7418M4
Efficiencies									
Power hi-lo	✓	✓	✓	✓	✓	✓	✓	✓	✓
Power backrest						✓	✓	✓	✓
Manual backrest	✓	✓	✓	✓	✓				
Foot control	✓	✓	✓	✓	✓	✓	✓	✓	✓
Front storage drawer	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reversible side drawer	✓	✓	✓			✓	✓		
Plastic drain pan	✓	✓	✓	✓	✓	✓	✓	✓	✓
Paper roll holder	✓	✓	✓	✓	✓	✓	✓	✓	✓
Foldaway stirrups		✓	✓	✓	✓	✓	✓	✓	✓
Pelvic tilt 7°			✓		✓		✓		✓
Infection Prevention									
OneTouch WheelBase® System				✓	✓			✓	✓
Classic seamless upholstery with PreFixx® Protective Finish	✓	✓	✓	✓	✓	✓	✓	✓	✓
Easy-clean powder coated enamel surfaces	✓	✓	✓	✓	✓	✓	✓	✓	✓
Patient Experience									
Drawer warmer			✓		✓		✓		✓
Accessibility									
500lb (226.8kg) patient weight capacity	✓	✓	✓	✓	✓	✓	✓	✓	✓
ADA Accessible 19" (48.3 cm)	✓	✓	✓	✓	✓	✓	✓	✓	✓
Electrical									
Hospital grade electrical receptacle			✓		✓		✓		✓
UL60601-1 and CAN/CSA-22.2 No. 601.1 classification	✓	✓	✓	✓	✓	✓	✓	✓	✓
Quality									
3-year limited warranty	✓	✓	✓	✓	✓	✓	✓	✓	✓
Durable galvanized steel construction	✓	✓	✓	✓	✓	✓	✓	✓	✓

Specifications

- Height: min. 19" (48.3 cm) at seat max. 32" (81.3 cm)
- Width: 28" (71.12 cm)
- Length: 58"-76" with foot section extended (147.3 cm – 193.0 cm)
- Leg rest: 18" W x 18" D (45.7 cm W x 45.7 cm D)
- Storage Drawer (Front): 18.0" W x 12.0" D x 3.0" H (46 cm x 30 cm x 8 cm)
- Storage Drawer (Side): 17.5" W x 21.0" D x 5.0" H (44 cm x 53 cm x 13 cm) (model dependent)
- Single column electro-magnetic actuator
- Electrical: 120 VAC, 60 HZ
- Amp range: 5 to 10.5 amp, depending on model
- Weight capacity: 500lb (226.8kg)
- Net weight: 325 lbs. (147kg) for 3000 series & 350 lbs. (159kg) for 3500 series
- Shipping weight: 365 lbs. (166kg) for 3000 series & 390 lbs. (177kg) for 3500 series
- Shipping dimensions: 30" W x 64" L x 31" H (76.2 cm x 162.6 cm x 78.7 cm)

Accessories

- 7419M1– Armrests
- 7419M2– Hand control for 3000 series
- 7419M3– Hand control for 3500 series

The OneTouch WheelBase[®] System

The OneTouch WheelBase[®] System makes moving our ProGlide[™] Series of products as easy as pushing a pedal. The integrated central-locking wheelbase with an intuitive foot pedal enables staff to move the equipment easily and quickly for cleaning, then lock it solidly into place for exams.

Model 7418M4
shown with optional
armrests



Upholstery & Colors

Every FusionONE and FusionONE+ model is available with Classic or Premium upholstery. Experience enhanced comfort, elegance, and softness in 16 stylish colors. The Premium upholstery model will elevate the aesthetics and comfort of any space.

Colors

							
Soft Linen	Smoky Cashmere	Warm Sand	Creamy Latte	Chocolate Truffle	Lemon Meringue	Cherry Blossom	Morning Fog
							
True Graphite	Classic Black	Fine Wine	Blue Skies	Lakeside Blue	Twilight Blue	Mint Leaf	Deep Forest