

Fusion**ONE**™

Power Hi-Lo, Manual Back Exam Table with Foot Control & Drain Pan

Model 7417M3

For examinations

With the accessibility of a power table and the affordability of a box table, the FusionONE Exam Table offers the features healthcare providers need most.

This hi-lo exam table enables simple ADA access and ease of examination for patients of all ages and abilities.



Key features at a glance

- › Power hi-lo and manual backrest
- › One-function foot control with hi-lo
- › Single column electro-magnetic actuator
- › 18"-21" paper roll holder (53.3cm)
- › 500lb (226.8kg) patient weight capacity
- › ADA Accessible 19" (48.3 cm)
- › Pelvic tilt (7°)
- › Drain pan: plastic 1.5" deep (3.81cm)
- › Large, reversible side drawer and large front drawer
- › Removable, seamless upholstered top with PreFixx® Protective Finish
- › Durable, galvanized steel construction
- › Hospital grade electrical receptacle with child safety cover
- › Heavy-gauge steel extension leg pull-out to accommodate most patient heights
- › Stirrups with four lateral positions, adjustable length
- › Easy-clean powder coat enamel surfaces
- › 3-year limited warranty

Specifications

- › Height: min. 19" (48.3 cm) at seat max. 32" (81.3 cm)
- › Width: 28" (71.12 cm)
- › Length: 58"-76" with foot section extended (147.3 cm – 193.0 cm)
- › Leg rest: 18" W x 18" D (45.7 cm W x 45.7 cm D)
- › Storage Drawer (Front) 18.0" W x 12.0" D x 3.0" H (46 cm x 30 cm x 8 cm)
- › Storage Drawer (Side) 17.5" W x 21.0" D x 5.0" H (44 cm x 53 cm x 13 cm)
- › Single column electro-magnetic actuator
- › Electrical: 120 VAC, 60 HZ
- › Weight capacity: 500lb (226.8kg)
- › Net weight: 300lb (136kg)
- › Shipping weight: 350lb (159kg)
- › Shipping dimensions: 30" W x 64" L x 31" H (76 cm x 163 cm x 79 cm)



Model 7417M3
shown with
options

Accessories

- › 7419M2 – Hand Control
- › 7419M1 – Armrests

Colors

			
Soft Linen	Smoky Cashmere	Warm Sand	Creamy Latte
			
True Graphite	Classic Black	Fine Wine	Blue Skies
			
Chocolate Truffle	Lemon Meringue	Cherry Blossom	Morning Fog
			
Lakeside Blue	Twilight Blue	Mint Leaf	Deep Forest



PRODUCTS ARE MADE IN THE USA

4767 Broadmoor Ave. SE, Suite #5, Kentwood, MI 49512

Phone: 616-514-1140 | Fax: 616-514-1141 | Toll Free: 877-514-1140 | www.medicus-health.com