

## FusionONE™ ProGlide™

Model 7417M4

### For examinations

With the accessibility of a power table and the affordability of a box table, the FusionONE ProGlide™ offers the features healthcare providers need most.

This hi-lo exam table with stirrups enables simple ADA access and ease of examination for patients of all ages and abilities.

The OneTouch WheelBase® System makes moving our ProGlide™ Series of products as easy as pushing a pedal. The integrated central-locking wheelbase with an intuitive foot pedal enables staff to move the equipment easily and quickly for cleaning, then lock it solidly into place for exams.



### Key features at a glance

- › OneTouch WheelBase® System
- › Power hi-lo and manual backrest
- › One-function foot control with hi-lo
- › Single column electro-magnetic actuator
- › 18"-21" paper roll holder (53.3cm)
- › 500lb (226.8kg) patient weight capacity
- › ADA Accessible 19" (48.3 cm)
- › Large front drawer
- › Drain pan: plastic 1.5' deep (3.81cm)
- › Removable, seamless upholstered top with PreFixx® Protective Finish
- › Heavy-gauge steel extension leg pull-out to accommodate most patient heights
- › Stirrups with four lateral positions, adjustable length
- › Durable, galvanized steel construction
- › Easy-clean powder coat enamel surfaces
- › 3-year limited warranty

## Specifications

- › Height: min. 19" (48.3 cm)  
at seat max. 32" (81.3 cm)
- › Width: 28" (71.12 cm)
- › Length: 58"-76" with foot section  
extended (147.3 cm – 193.0 cm)
- › Leg rest: 18" W x 18" D (45.7 cm W  
x 45.7 cm D)
- › Storage Drawer (Front): 18.0" W x 12.0"  
D x 3.0" H (46 cm x 30 cm x 8 cm)
- › Single column electro-magnetic actuator
- › Electrical: 120VAC, 60 HZ
- › Weight capacity: 500lb (226.8kg)
- › Net weight: 300lb (136kg)
- › Shipping weight: 350lb (159kg)
- › Shipping dimensions: 30" W x  
64" L x 31" H (76 cm x 163 cm x 79 cm)

## Accessories

- › 7419M2 – Hand Control
- › 7419M1 – Armrests



Model **7417M4**  
shown with options

## Upholstery and Colors

Every model is available with Classic  
or Premium upholstery and in 16 colors.

			
Soft Linen	Smoky Cashmere	Warm Sand	Creamy Latte
			
True Graphite	Classic Black	Fine Wine	Blue Skies
			
Chocolate Truffle	Lemon Meringue	Cherry Blossom	Morning Fog
			
Lakeside Blue	Twilight Blue	Mint Leaf	Deep Forest



PRODUCTS ARE MADE IN THE USA

4767 Broadmoor Ave. SE, Suite #5, Kentwood, MI 49512

Phone: 616-514-1140 | Fax: 616-514-1141 | Toll Free: 877-514-1140 | [www.medicus-health.com](http://www.medicus-health.com)