

FusionONE+™

Power Hi-Lo, Power Back Exam Table with Foot Control, Drain Pan & Stirrups

For examinations

Fusing the accessibility of a power table and the affordability of a fixed-height box table, the FusionONE+ Power Exam Table offers the features healthcare providers need most.

Featuring a power back, this hi-lo exam table can accommodate multiple specialties. It enables simple ADA access and ease of examination for patients of all ages and abilities. The power backrest makes it easy to lie patients back and sit them up safely and comfortably.



Model 7418M2

PROTECTS

- › patients from infection
- › facilities from damage
- › staff from injury

Key features at a glance

- › Power hi-lo and power backrest
- › Two-function foot control
- › 18"-21" (45.7 cm-53.3 cm) paper roll holder
- › 500 lb. (226.8 kg) patient weight capacity
- › ADA Accessible 19" (48.3 cm)
- › Pelvic tilt (7°)
- › Drain pan: plastic 1.5" deep (3.8 cm)
- › Large, reversible side drawer and large front drawer
- › Removable, seamless upholstered top with PreFixx® Protective Finish
- › Durable, galvanized steel construction
- › Electrical receptacle
- › Drawer warmer
- › Heavy-gauge steel pull-out leg extension
- › Foldaway stirrups with four lateral positions and adjustable length
- › Easy-clean powder coat enamel surfaces
- › 3-year limited warranty

Specifications

- › Height: min. 19" (48.3 cm) at seat
max. 32" (81.3 cm)
- › Width: 28" (71.1 cm)
- › Length: 58"-76" (147.3 cm-193 cm)
with leg rest extended
- › Leg rest: 18" W x 18" D (45.7 cm W x 45.7 cm D)
- › Storage Drawer (Front): 18" W x 12" D x 3" H
(45.7 cm x 30.5 cm x 7.6 cm)
- › Storage Drawer (Side): 17.5" W x 21" D x 5" H
(44.5 cm x 53.3 cm x 12.7 cm)
- › Single-column, electro-magnetic actuator
- › Electrical: 120 VAC 10.5 AMP, 60Hz
- › Weight capacity: 500 lbs. (226.8kg)
- › Net weight: 350 lbs. (159kg)
- › Shipping weight: 390 lbs. (177kg)
- › Shipping dimensions: 30" W x 64" L x 31" H
(76.2 cm x 162.6 cm x 78.7 cm)

Accessories

- › **7419M3** – Hand Control
- › **7419M1** – Armrests

Colors



PRODUCTS ARE MADE IN THE USA