

## Fusion**ONE+**<sup>TM</sup> ProGlide<sup>TM</sup>

Model 7418M4

### For examinations

Fusing the accessibility of a power table and the affordability of a fixed-height box table, the FusionONE+ ProGlide offers the features healthcare providers need most.

Featuring a power back, this hi-lo exam table can accommodate multiple specialties. It enables simple ADA access and ease of examination for patients of all ages and abilities. The power backrest makes it easy to lie patients back and sit them up safely and comfortably.

The OneTouch WheelBase® System makes moving our ProGlide™ Series of products as easy as pushing a pedal. The integrated central-locking wheelbase with an intuitive foot pedal enables staff to move the equipment easily and quickly for cleaning, then lock it solidly into place for exams.



### Key features at a glance

- › OneTouch WheelBase® System
- › Power hi-lo and power backrest
- › Two-function foot control
- › 18"-21" (45.7 cm-53.3 cm) paper roll holder
- › 500 lb. (226.8 kg) patient weight capacity
- › ADA Accessible 19" (48.3 cm)
- › Pelvic tilt (7°)
- › Drain pan: plastic 1.5" deep (3.8 cm)
- › Large front drawer
- › Removable, seamless upholstered top with PreFixx® Protective Finish
- › Durable, galvanized steel construction
- › Electrical receptacle
- › Heavy-gauge steel pull-out leg extension
- › Foldaway stirrups with four lateral positions and adjustable length
- › Easy-clean powder coat enamel surfaces
- › 3-year limited warranty

## Specifications

- › Height: min. 19" (48.3 cm)  
at seat max. 32" (81.3 cm)
- › Width: 28" (71.1 cm)
- › Length: 58"-76" (147.3 cm-193 cm) with  
leg rest extended
- › Leg rest: 18" W x 18" D (45.7 cm W  
x 45.7 cm D)
- › Storage Drawer (Front): 18" W x 12" D x  
3" H (45.7 cm x 30.5 cm x 7.6 cm)
- › Single-column, electro-magnetic actuator
- › Electrical: 120 VAC 10.5 AMP, 60Hz
- › Weight capacity: 500 lbs. (226.8kg)
- › Net weight: 350 lbs. (159kg)
- › Shipping weight: 390 lbs. (177kg)
- › Shipping dimensions: 30" W x 64" L x 31" H  
(76.2 cm x 162.6 cm x 78.7 cm)



Model **7418M4**  
shown with options

## Accessories

- › 7419M3 – Hand Control
- › 7419M1 – Armrests

## Colors

			
Soft Linen	Smoky Cashmere	Warm Sand	Creamy Latte
			
True Graphite	Classic Black	Fine Wine	Blue Skies
			
Chocolate Truffle	Lemon Meringue	Cherry Blossom	Morning Fog
			
Lakeside Blue	Twilight Blue	Mint Leaf	Deep Forest



PRODUCTS ARE MADE IN THE USA

4767 Broadmoor Ave. SE, Suite #5, Kentwood, MI 49512

Phone: 616-514-1140 | Fax: 616-514-1141 | Toll Free: 877-514-1140 | [www.medicus-health.com](http://www.medicus-health.com)