

proglide300 power exam table

Power Hi-Lo, Manual Back Exam Table with OneTouch WheelBase®, OneTouch Patient Positioning® System, Hand and Foot Control

For examinations

The proglide300 is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers. Designed for exams, the proglide300 meets ADA and U.S. Access Board regulations and accommodates the needs of most patients.



Model 7420M3

PROTECTS

- patients from infection
- facilities from damage
- staff from injury

Key features at a glance

- › OneTouch WheelBase® integrated central-locking caster system
- › Hand control with OneTouch Patient Positioning®, 2 programmable positions & return-to-home
- › One function foot control
- › Power hi-lo with smooth moving articulating back adjustment
- › Dual electro-mechanical actuator system for efficient speed
- › 18"-21" paper roll holder (53.3cm)
- › 650lb (294.84kg) patient weight capacity
- › Low-access ADA height 18.5" (±0.5") - 46.9cm (±1.27cm)
- › Pelvic tilt (7°)
- › Drain pan: stainless steel 1.5" deep (3.81cm)
- › Front storage drawer
- › Side storage drawer – reversible
- › Removable, seamless upholstery with Prefixx® Protective Finish
- › Stirrups with four lateral positions, adjustable length
- › Front leg cushion for patient comfort
- › Heavy gauge steel extension leg pull-out to accommodate all patient heights
- › Drawer warmer
- › Hospital grade receptacle 110V with child safety cover (patient left)
- › 3 year limited warranty
- › Durable galvanized steel construction
- › Easy-clean powder coat enamel surfaces



Specifications

- › Height: min. 18.5" ($\pm .5$ ") 46.9cm (± 1.27 cm) at seat max. 37" at seat (93.9cm)
- › Width: 28" (71.12cm)
- › Length: 57" – 76" with foot section extended (144.7cm - 193cm)
- › Leg rest: 19" W x 13" D (48.26cm W x 33cm D)
- › Front storage drawer: 19" W x 12" D (48.26cm W x 30.48cm D)
- › Stirrups: four lateral positions, adjustable length
- › Dual electro-mechanical actuator system for efficient speed
- › Electrical: 120 VAC, 60 HZ
- › Net weight: 470lb (213kg)
- › Shipping weight: 510lb (231kg)
- › Shipping dimensions: L64xW30xH31 (163cm x 76cm x 79cm)

Accessories

- › 7422M1 – Patient assist arms complete (pair) field installable
- › 7422M2 – Universal Accessory Rail-Seat Section
- › 7422M3 – Removable stainless steel drain pan W/spout
- › 7422M4 – 110v drawer warmer (factory installed)

Colors



4767 Broadmoor Ave. SE, Suite #5, Kentwood, MI 49512

Phone: 616-514-1140 | Fax: 616-514-1141 | Toll Free: 877-514-1140 | www.medicus-health.com