

power100+ power exam table

Power Hi-Lo, Power Back Exam Table with Foot Control

Model 7421M1

For examinations

The power100+ is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers.

Designed for exams, the power100+ meets ADA and U.S. Access Board regulations and accommodates the needs of most patients.



PROTECTS > patients from infection > facilities from damage > staff from injury

Key features at a glance

- > Power hi-lo + power backrest
- > Two-function foot control with hi-lo and power back
- > Dual electro-mechanical actuator system for efficient speed
- > 18"-21" paper roll holder (53.3cm)
- > 650lb (294.84kg) patient weight capacity
- > Low-access ADA height 18.5" (±0.5") - 46.9cm (±1.27cm)
- > Pelvic tilt (7°)
- > Drain pan: stainless steel 1.5" deep (3.81cm)
- > Large front storage drawer
- > Removable, seamless upholstery with Prefixx® Protective Finish
- > Stirrups with four lateral positions, adjustable length
- > Front leg cushion for patient comfort
- > Heavy gauge steel extension leg pull-out to accommodate all patient heights
- > 3 year limited warranty
- > Durable galvanized steel construction
- > Easy-clean powder coat enamel surfaces

Wireless Hand & Foot
Controls Available



Specifications

- › Height: min. 18.5" (± .5") 46.9cm (±1.27cm)
at seat max. 37" at seat (93.9cm)
- › Width: 28" (71.12cm)
- › Length: 57" – 76" with foot section
extended (144.7cm - 193cm)
- › Leg rest: 19" W x 13" D
(48.26cm W x 33cm D)
- › Front storage drawer: 19" W x 12" D
(48.26cm W x 30.48cm D)
- › Stirrups: four lateral positions, adjustable
length
- › Dual electro-mechanical actuator system
for efficient speed
- › Electrical: 120 VAC, 60 HZ
- › Net weight: 405lb (184kg)
- › Shipping weight: 445lb (201kg)
- › Shipping dimensions: L64 x W30 x H31
(163cm x 76cm x 79cm)

Accessories

- › **7422M1** – Patient assist arms complete (pair) field installable
- › **7422M2** – Universal Accessory Rail-Seat Section
- › **7422M3** – Removable stainless steel drain pan W/spout
- › **7422M4** – 110v drawer warmer (factory installed)
- › **7422M5** – Two-function wireless foot control
- › **7422M6** – Two-function wireless hand control

Colors



PRODUCTS ARE MADE IN THE USA