

proglide300⁺ power exam table

Power Hi-Lo, Power Back Exam Table with OneTouch WheelBase[®] System, OneTouch Patient Positioning[®] System, Hand and Foot Control

Model 7421M3



For examinations

The proglide300+ is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers. Designed for exams, the proglide300+ meets ADA and U.S. Access Board regulations and accommodates the needs of most patients.

PROTECTS

➤ patients from infection ➤ facilities from damage ➤ staff from injury

Key features at a glance

- › OneTouch WheelBase[®] integrated central-locking caster system
- › Power hi-lo and power backrest
- › Two function foot control with hi-lo and power back
- › Hand control with OneTouch Patient Positioning[®], 2 programmable positions & return-to-home
- › Dual electro-mechanical actuator system for efficient speed
- › 18"-21" paper roll holder (53.3cm)
- › 650lb (294.84kg) patient weight capacity
- › Low-access ADA height 18.5" (±0.5") - 46.9cm (±1.27cm)
- › Pelvic tilt (7°)
- › Drain pan: stainless steel 1.5" deep (3.81cm)
- › Front storage drawer
- › Side storage drawer – reversible
- › Removable, seamless upholstery with Prefixx[®] Protective Finish
- › Stirrups with four lateral positions, adjustable length
- › Front leg cushion for patient comfort
- › Heavy gauge steel extension leg pull-out to accommodate all patient heights
- › Drawer warmer
- › Hospital grade receptacle 110V with child safety cover (patient left)
- › 3 year limited warranty
- › Durable galvanized steel construction
- › Easy-clean powder coat enamel surfaces

Wireless Hand & Foot
Controls Available





Specifications

- › Height: min. 18.5" ($\pm .5$ ") 46.9cm (± 1.27 cm)
at seat max. 37" at seat (93.9cm)
- › Width: 28" (71.12cm)
- › Length: 57" – 76" with foot section
extended (144.7cm - 193cm)
- › Leg rest: 19" W x 13" D
(48.26cm W x 33cm D)
- › Front storage drawer: 19" W x 12" D
(48.26cm W x 30.48cm D)
- › Stirrups: four lateral positions,
adjustable length
- › Dual electro-mechanical actuator
system for efficient speed
- › Electrical: 120 VAC, 60 HZ
- › International model 4040-650-301:
220/240V, 50/60Hz
- › Net weight: 445lb (202kg)
- › Shipping weight: 485lb (220kg)
- › Shipping dimensions: L64 x W30 x H31
(163cm x 76cm x 79cm)

Accessories

- › 7422M1 – Patient assist arms complete (pair) field installable
- › 7422M2 – Universal Accessory Rail-Seat Section
- › 7422M3 – Removable stainless steel drain pan W/spout
- › 7422M4 – 110v drawer warmer (factory installed)
- › 7422M5 – Two-function wireless foot control
- › 7422M6 – Two-function wireless hand control

Colors



PRODUCTS ARE MADE IN THE USA

4767 Broadmoor Ave. SE, Suite #5, Kentwood, MI 49512

Phone: 616-514-1140 | Fax: 616-514-1141 | Toll Free: 877-514-1140 | www.medicus-health.com