

power4010

head-centric procedure chair

Procedure Chair with OneTouch Patient Positioning® System

For head-centric practices

The power4010 is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers, such as ENT, Dermatology, and Plastic Surgery professionals. Designed for head-centric procedures, the power4010 meets ADA and U.S. Access Board regulations and accommodates the needs of most patients.



Model 7425M1

PROTECTS

- > patients from infection
- > facilities from damage
- > staff from injury

Key features at a glance

- > Chair designed for head-centric practices
- > Slim back cushion for accessibility
- > Double-articulating headrest with height adjustability and interchangeable headrest options
- > 650 lb. (294.84 kg) patient weight capacity
- > Low access ADA height 18.5" (±0.5") 46.9cm (±1.27cm)
- > 4 power functions: hi-lo, back, tilt, leg lift & extension
- > Hand control with OneTouch Patient Positioning® System, 2 programmable positions and auto-return
- > Power seat that tilts up
- > Power leg lift for patient comfort
- > Concealed 18" paper roll holder
- > Trendelenberg, Lithotomy & Semi-Fowler positions
- > Removable, seamless upholstery with Prefixx® Protective Finish
- > Hospital-grade electrical receptacle with child safety cover (patient left)

Hand
Control



Specifications

- › Height: min. 18.5" (±.5") 46.9cm (±1.27cm) at seat max. 37" at seat 93.9cm
- › Seat width: 28" (71.12cm)
- › Back width: 24.5" (62.23cm)
- › Length: 70" – 80" with foot section extended (177.8cm – 203.2cm)
- › Leg rest: 28" W x 13" D (71.12cm W x 33.02cm D)
- › Dual electro-mechanical actuator system for efficient speed
- › Electrical: 120 VAC, 60HZ, 5 AMP
- › Net weight: 410lb (185kg)
- › Shipping weight: 450lb (204kg)
- › Shipping dimensions: L64 x W30 x H41 (163cm x 76cm x 104cm)

Accessories

- › 7428M1 – Adjustable / self-leveling arms (pair)
- › 7428M3 – Removable U-shaped headrest
- › 7428M4 – Removable rectangular headrest



Shown with options (7428M1)

Colors

			
Soft Linen	Smoky Cashmere	Warm Sand	Creamy Latte
			
Chocolate Truffle	Lemon Meringue	Cherry Blossom	Morning Fog
			
True Graphite	Classic Black	Fine Wine	Blue Skies
			
Lakeside Blue	Twilight Blue	Mint Leaf	Deep Forest



PRODUCTS ARE MADE IN THE USA