

## proglide4010

head centric procedure chair

### Procedure Chair with OneTouch WheelBase® System, OneTouch Patient Positioning® System, Hand and Foot Control

#### For head-centric practices

The proglide4010 is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers, such as ENT, Dermatology, and Plastic Surgery professionals.

Designed for head-centric procedures, the proglide4010 meets ADA and U.S. Access Board regulations and accommodates the needs of most patients.



#### PROTECTS

- > patients from infection
- > facilities from damage
- > staff from injury

#### Key features at a glance

- > Chair designed for head-centric practices
- > OneTouch Wheelbase®
- > 650 lb. (294.84 kg) patient weight capacity
- > Low access ADA height 18.5" (±0.5") 46.9cm (±1.27cm)
- > 4 Power functions: hi-lo, back, tilt, leg lift & extension
- > Hand control with OneTouch Patient Positioning® System, 2 programmable positions & return-to-home
- > Foot control with OneTouch Patient Positioning® System, 2 programmable positions and return-to-home
- > Double-articulating headrest with height adjustability and interchangeable headrest options
- > Slim back cushion for accessibility
- > Removable, seamless upholstery top with Prefixx® Protective Finish
- > Trendelenberg, Lithotomy & Semi-Fowler positions
- > Hospital-Grade electrical receptacle with child safety cover (patient left)
- > Concealed 18" paper roll holder
- > 3 year limited warranty

#### Hand Control



## Specifications

- › Height: min. 18.5" (± .5") 46.9cm (±1.27cm)  
at seat max. 37" at seat (93.9cm)
- › Seat width: 28" (71.12cm)
- › Back width: 24.5" (62.23cm)
- › Length: 70" – 80" with foot section  
extended (177.8cm – 203.2cm)
- › Leg rest: 28" W x 13" D (71.12cm  
W x 33.02cm D)
- › Dual electro-mechanical actuator system  
for efficient speed
- › Electrical: 120 VAC, 60HZ , 5 AMP
- › Net weight: 450lb (204kg)
- › Shipping weight: 490lb (222kg)
- › Shipping dimensions: L64 x W30 x H41  
(163cm x 76cm x 104cm)

## Accessories

- › 7428M1 – Adjustable / self-leveling arms (pair)
- › 7428M3 – Removable U-shaped headrest
- › 7428M4 – Removable rectangular headrest



Shown with  
options (7428M1)

## Colors

			
Soft Linen	Smoky Cashmere	Warm Sand	Creamy Latte
			
Chocolate Truffle	Lemon Meringue	Cherry Blossom	Morning Fog
			
True Graphite	Classic Black	Fine Wine	Blue Skies
			
Lakeside Blue	Twilight Blue	Mint Leaf	Deep Forest



PRODUCTS ARE MADE IN THE USA