

proglide4011

ultra procedure chair

Ultra Procedure Chair with OneTouch WheelBase®, OneTouch Patient Positioning® System Hand and Foot Control

For procedures

The proglide4011 is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers, such as OB/GYN and Oncology professionals.

Designed for procedures, the proglide4011 meets ADA and U.S. Access Board regulations and accommodates the needs of most patients.



Model 7426M3

PROTECTS

- patients from infection
- facilities from damage
- staff from injury

Key features at a glance

- › OneTouch WheelBase® integrated central-locking caster system
- › 650 lb. (294.84kg) patient weight capacity
- › Low access ADA height 18.5" (±0.5") 46.9cm (±1.27cm)
- › 4 power functions: hi-lo, back, tilt, leg lift & extension
- › Hand control with OneTouch Patient Positioning® System, 2 programmable positions & return-to-home
- › Foot control with OneTouch Patient Positioning® System, 2 programmable positions & return-to-home
- › Removable, seamless upholstered top with Prefix® Protective Finish
- › Double-articulating headrest with adjustable height
- › Stirrups
- › Interchangeable headrest options
- › Stainless steel drain pan
- › Trendelenburg, Lithotomy, Semi-Fowler positions
- › Hospital-grade electrical receptacle with child safety cover (patient left)
- › Concealed paper roll holder
- › 3 year limited warranty

Hand Control



Specifications

- › Height: min. 18.5" (± .5") 46.9cm (±1.27cm)
at seat max. 37" at seat (93.9cm)
- › Seat width: 28" (71.12cm)
- › Back width: 26.5" (67.31cm)
- › Length: 73"-80" with foot section
extended (185.42cm – 203.2cm)
- › Leg rest: 28" W x 11.5" D
(71.12cm W x 29.21cm D)
- › Stirrups: four lateral
positions, adjustable length
- › Dual electro-mechanical actuator system
for efficient speed
- › Electrical: 120 VAC, 60HZ
- › Net weight: 495lb (224kg)
- › Shipping weight: 535lb (242kg)
- › Shipping dimensions: L64 x W30 x H41
(163cm x 76cm x 104cm)

Accessories

- › 7428M2 – Adjustable / self-leveling arms (pair)
- › 7428M3 – Removable U-shaped headrest
- › 7428M4 – Removable rectangular headrest



Shown with
options (7428M2)

Colors

Soft Linen	Smoky Cashmere	Warm Sand	Creamy Latte
Chocolate Truffle	Lemon Meringue	Cherry Blossom	Morning Fog
True Graphite	Classic Black	Fine Wine	Blue Skies
Lakeside Blue	Twilight Blue	Mint Leaf	Deep Forest



PRODUCTS ARE MADE IN THE USA