

7429M

FEATURES & ANATOMY



SPECIFICATIONS

Maximum Patient Weight	350 lb. (159 kg.)
Seat Height, Low	25.1 in. (64 cm.)
Seat Height, High	33 in. (83.82 cm.)
Seat Cushion Height, Front (tapered)	3.25 in. (8.26 cm.)
Seat Cushion Height, Back (tapered)	2.25 in. (5.08 cm.)
Back Cushion Height	33 in. (83.82 cm.)
Seat Width	24 in. (61 cm.)
Back Width (no extension)	15 in. (38 cm.)
Back Width (w/ optional extension)	24 in. (61 cm.)
Radiolucent Window (width x height)	13 in. x 26 in. (33 cm. x 65 cm.)
Overall Width (rails up)	30.5 in. (77.47 cm.)
Patient Surface Length (standard)	77 in. (195.58 cm.)
Patient Surface Length (w/ manual folding footrest)	70.4 in. (179 cm.)
Overall Chair Height, Low	55.5 in. (141 cm.)