

7430M

FEATURES & ANATOMY



SPECIFICATIONS

	7430M1	7430M2	7430M3	7430M4	7430M5
Maximum Patient Weight	500 lb. (227 kg)	500 lb. (227 kg)	500 lb. (227 kg)	500 lb. (227 kg)	500 lb. (227 kg)
Seat Height, Low	24.3" (61.7 cm)	24.3" (61.7 cm)	24.3" (61.7 cm)	24.3" (61.7 cm)	24.4" (62 cm)
Seat Height, High	32.1" (81.5 cm)	32.1" (81.5 cm)	32.1" (81.5 cm)	32.1" (81.5 cm)	40.1" (101.9 cm)
Back Height	33.75" (85.73 cm)	33.75" (85.73 cm)	33.75" (85.73 cm)	33.75" (85.73 cm)	33.75" (85.73 cm)
Side Rail Height, Up (seat pan to rail top)	13.5" (34.3 cm)	13.5" (34.3 cm)	13.5" (34.3 cm)	13.5" (34.3 cm)	13.5" (34.3 cm)
Overall Chair Height, Low	54" (137 cm)	54" (137 cm)	54" (137 cm)	54" (137 cm)	53" (135 cm)
Overall Chair Height, High	62" (157 cm)	62" (157 cm)	62" (157 cm)	62" (157 cm)	70" (178 cm)
Seat Width	24" (61 cm)	24" (61 cm)	28" (71 cm)	28" (71 cm)	28" (71 cm)
Overall Width (rails up)	30" (76 cm)	30" (76 cm)	34" (86 cm)	34" (86 cm)	34" (86.4 cm)
Overall Width (rails down)	24" (61 cm)	24" (61 cm)	28" (71 cm)	28" (71 cm)	28" (71 cm)
Patient Surface Length (without "F" option)	77" (195 cm)	77" (195 cm)	77" (195 cm)	77" (195 cm)	77" (195 cm)
Side Rail Length	27.6" (70.1 cm)	27.6" (70.1 cm)	27.6" (70.1 cm)	27.6" (70.1 cm)	27.6" (70.1 cm)
Backrest Articulation Range	0° - 75°	0° - 75°	0° - 75°	0° - 75°	0° - 75°
Trendelenburg	5°	5°	5°	5°	5°
Radiolucent Window (width x height)	20" X 25.5" (51 cm X 65 cm)	20" X 25.5" (51 cm X 65 cm)	20" X 25.5" (51 cm X 65 cm)	20" X 25.5" (51 cm X 65 cm)	20" X 25.5" (51 cm X 65 cm)
Wheelbase (length X width)	30" x 20" (76 cm x 51 cm)	30" x 20" (76 cm x 51 cm)	30" x 20" (76 cm x 51 cm)	30" x 20" (76 cm x 51 cm)	30" x 20" (76 cm x 51 cm)