




**Traveling Arms Recliner**  
Item# 2502M6

#### SPECIFICATIONS:

- Traveling Arms
- 3-Position Recline: Up-right, Semi, and Full
- 6" Adjustable Arm Height
- Extra-wide Padded "L" Arms



Seat Height	Patient Width	Weight Capacity	Overall Dimensions
25"	25"	325 lbs	39"W x 42"D (Upright)

#### READY TO SHIP PRODUCT: 12RMAW

This chair is currently not available as a Ready To Ship item.

#### FEATURES

This unique X-Wide reclining chair allows the arms to travel back with the patient as they are reclined. The arms are fully adjustable up/down as well as in/out.

#### ACCESSORIES



Matching Head Cover



Flip-up Side Trays

#### FINISH SELECTIONS

To view our **Vinyl Color Selection**, please visit:  
<https://www.medicus-health.com/wheeled-reclining-treatment-chair-extra-wide.html>

#### WARRANTY INFORMATION

To view our **Limited Warranty Policy**, please visit: <https://www.medicus-health.com/wheeled-reclining-treatment-chair-extra-wide.html>

#### PARTS INFORMATION

To view our **Parts List** for this product, please visit:  
[www.medicus-health.com/wheeled-reclining-treatment-chair-extra-wide.html](https://www.medicus-health.com/wheeled-reclining-treatment-chair-extra-wide.html)

#### CERTIFICATION(S)

As an option we offer CAL-133 protection, a flammability test procedure for seating furniture for use in Public Occupancies that has been adopted by many states. The test results show that the product passed both Criteria Q and Criteria B for the California TB-133 Burn Test.

