

Extended hand protection for a variety of jobs

Thick design: gloves provide reliable protection with its thick composition

Extended cuff: These nitrile gloves also offers an expanded shield over the wrist and forearm

Food handling compliance: Suitable for a variety of food handling applications



Industries

- Life Sciences

Recommended For

- Maintaining lab furniture and equipment
- Sample taking and processing

TECHNICAL DATA SHEET

	Product Information
External Glove Surface	Textured Fingers
Packaging Overview	50 gloves per dispenser 10 dispensers per case 500 gloves per case
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Country Of Origin	Malaysia
Antistatic	Not Tested

PHYSICAL PROPERTIES

	Typical Values		Testing Method
Length (mm/inches)	300 / 11.8		ASTM D3767,EN 420
Freedom from Holes (Inspection level I)	1.5 AQL		ASTM D5151,EN 455-1
Palm Thickness (mm/mils)	0.15/5.9		ASTM D3767,EN 420
Finger Thickness (mm/mils)	0.19/7.5		ASTM D3767,EN 420
	BEFORE AGING	AFTER AGING	
Ultimate Tensile Strength (MPa)	27	14	ASTM D412 & D573
Elongation at Break (%)	≥ 500	≥ 400	ASTM D412

ORDERING INFORMATION

Size	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	7456M1	7456M2	7456M3	7456M4

Performance Standards and Regulatory Compliance

