

Black Dragon® Zero

Black Nitrile Exam Glove with Textured Fingertips

Tough, durable protection in a black color to hide dirt, grease and stains

Black nitrile: Black Dragon Zero are black nitrile medical exam gloves that mask stains, grease, and dirt and provides visual contrast

Secure grip: With textured fingertips, the gloves help provide a secure grip

Elevated comfort and strength: These gloves are designed with a soft, durable nitrile formulation



Industries

- Automotive
- Automotive Aftermarket
- Healthcare

Recommended For

- Inspection, selecting, checking parts
- Inspection of parts, equipment
- Picking and fastening components
- Equipment repair and maintenance
- General Purpose Auto Aftermarket
- Use of cleaning tools and robotic equipment
- Tire rotation & changes
- Standard, moderate risk examination procedures
- Protection from Type I latex allergy in HCW's or patients

TECHNICAL DATA SHEET

	Product Information
External Glove Surface	Textured Fingers
Audit Standards	EN ISO 13485:2012
Packaging Overview	100 gloves per dispenser 10 dispensers per case 1000 gloves per case
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Country Of Origin	Thailand
Antistatic	Not Tested

PHYSICAL PROPERTIES

	Typical Values		Testing Method
Length (mm/inches)	245 / 9.6		ASTM D3767,EN 420
Freedom from Holes (Inspection level I)	1.5 AQL		ASTM D5151,EN 455-1
Palm Thickness (mm/mils)	0.09 / 3.5		ASTM D3767,EN 420
Finger Thickness (mm/mils)	0.14 / 5.5		ASTM D3767,EN 420
	BEFORE AGING	AFTER AGING	
Ultimate Tensile Strength (MPa)	≥ 33	≥ 29	ASTM D412 & D573
Elongation at Break (%)	≥ 500	≥ 400	ASTM D412

ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)	XXL (10.5 - 11)
Product Code	7460M1	7460M2	7460M3	7460M4	7460M5	7460M6