

# ComfortGrip®

## Latex Exam Glove

# Affordable, comfortable hand protection with a secure grip

**Comfortable fit and secure grip:**

ComfortGrip® are comfortable, powder-free gloves with a fully textured surface to allow wearers in a range of industries to securely grip tools and materials

**Medical exam grade:** These disposable latex gloves can be used by healthcare and dental professionals in certain applications as they are medical exam grade certified

**Food handling compliance:** These gloves are also suitable for a wide range of food-handling applications

**Great value:** Hand protection you can trust at an excellent price point

**CAUTION:** This product contains natural rubber latex, which may cause allergic reactions

**Industries**

- Food Processing
- Life Sciences
- Healthcare

**Recommended For**

- All food processing, applications that have food contact
- Sample taking and processing
- Extra comfort during long, sustained procedures
- Standard patient examinations

TECHNICAL DATA SHEET

|                        | Product Information                                                                                      |
|------------------------|----------------------------------------------------------------------------------------------------------|
| External Glove Surface | Fully Textured                                                                                           |
| Audit Standards        | ISO 13485:2003                                                                                           |
| Packaging Overview     | 100 gloves per dispenser<br>10 dispensers per case<br>1000 gloves per case                               |
| Storage Instructions   | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition. |
| Country Of Origin      | Malaysia                                                                                                 |
| Product Segmentation   | Multi-Purpose                                                                                            |
| Antistatic             | Not Tested                                                                                               |

PHYSICAL PROPERTIES

|                                         | Typical Values |             | Testing Method      |
|-----------------------------------------|----------------|-------------|---------------------|
| Length (mm/inches)                      | 245 / 9.6      |             | ASTM D3767,EN 420   |
| Freedom from Holes (Inspection level I) | 1.5 AQL        |             | ASTM D5151,EN 455-1 |
| Palm Thickness (mm/mils)                | 0.13 / 5.1     |             | ASTM D3767,EN 420   |
| Finger Thickness (mm/mils)              | 0.16 / 6.3     |             | ASTM D3767,EN 420   |
|                                         | BEFORE AGING   | AFTER AGING |                     |
| Ultimate Tensile Strength (MPa)         | 22             | 20          | ASTM D412 & D573    |
| Elongation at Break (%)                 | ≥ 700          | ≥ 600       | ASTM D412           |

ORDERING INFORMATION

| Size         | XS (5.5 - 6) | S (6.5 - 7) | M (7.5 - 8) | L (8.5 - 9) | XL (9.5 - 10) |
|--------------|--------------|-------------|-------------|-------------|---------------|
| Product Code | 7465M1       | 7465M2      | 7465M3      | 7465M4      | 7465M5        |