

Diamond Grip Plus®

Latex Exam Glove

Latex exam gloves with exceptional fit and feel for extended wear

Secure grip: Diamond Grip Plus latex medical exam gloves are fully textured, for a secure grip in both wet or dry environments

Tactile sensitivity: Excellent tactile sensitivity ensures precise handling during patient care or when handling delicate products or laboratory instruments

Comfortable and elastic fit: Superior elasticity means gloves move with you for greater comfort

Outstanding quality: Exceeds standards for tensile strength, elasticity/elongation and other quality standards for average latex gloves

CAUTION: This product contains natural rubber latex, which may cause allergic reactions



Industries

- Food Processing
- Life Sciences
- Healthcare
- Utilities
- Warehousing

Recommended For

- Equipment, maintenance and repair
- All food processing, applications that have food contact
- Sample taking and processing
- Extra comfort during long, sustained procedures
- Standard patient examinations

TECHNICAL DATA SHEET

	Product Information
External Glove Surface	Fully Textured
Audit Standards	EN ISO 13485:2012
Packaging Overview	100 gloves per dispenser 10 dispensers per case 1000 gloves per case
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Country Of Origin	Malaysia
Antistatic	Not Tested

PHYSICAL PROPERTIES

	Typical Values		Testing Method
Length (mm/inches)	245 / 9.6		ASTM D3767,EN 420
Freedom from Holes (Inspection level I)	1.5 AQL		ASTM D5151,EN 455-1
Palm Thickness (mm/mils)	0.14 / 5.5		ASTM D3767,EN 420
Finger Thickness (mm/mils)	0.16 / 6.3		ASTM D3767,EN 420
	BEFORE AGING	AFTER AGING	
Ultimate Tensile Strength (MPa)	22	20	ASTM D412 & D573
Elongation at Break (%)	≥ 700	≥ 600	ASTM D412
Force at break (N)	9	9	EN 455-2

ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	7467M1	7467M2	7467M3	7467M4	7467M5