

# Powered Bariatric Recliner

## Comfort - Functionality - Durability

The Powered Bariatric Recliner is the perfect solution for larger patients who require a wider recliner with an increased weight capacity. This recliner, with a heavy-gauge steel frame, has a maximum weight capacity of 700 lb when

evenly distributed! A four-function hand pendant activates the motor, allowing caregivers to safely recline and reposition patients.



Multiple Positions for Enhanced Patient Comfort

### Features:

- Quick-release trigger mechanism on push handle allows recliner to swiftly move into Trendelenburg position
- Equipped with rechargeable battery pack for operation if the power goes out
- Four swivel casters facilitate transport and maneuverability in confined areas
- Durable, high impact side tables sit flush with chair arm to conveniently hold prep materials. Side table is smooth and features tapered edges for easy cleaning
- Legrest and backrest raise and lower independently, giving the caregiver complete control over positioning
- Comes complete with head pillow for extra support and patient comfort

### Dimensions

Overall Height .....	55" / 140 cm
Overall Width .....	38" / 97 cm
Overall Depth .....	59" / 150 cm
Width Between Arms .....	25.25" / 64 cm
Arm Height From Seat .....	8" / 20 cm
Seat Height .....	25.5" / 65 cm
Seat Width .....	25.75" / 65 cm
Seat Depth .....	24" / 61 cm
Depth Fully Reclined In Trendelenburg .....	71" / 180 cm
Width With Both Side Tables in Upright Position .....	58" / 147 cm
Maximum Weight Capacity Evenly Distributed .....	700 lb / 317.5 kg



### Available in Four Stocking Colors:

- FR588W427** Blue Ridge
- FR588W857** Jade
- FR588W863** Rosewood
- FR588W851** Doe Skin

(Meets California Technical Bulletin 133 Flammability Standards)



Blue Ridge (427)



Jade (857)



Rosewood (863)



Doe Skin (851)