

# 2959M1 - Deluxe Preferred Care Geriatric Recliner, Wide

The Preferred Care Recliner-Wide is favored for its 24 inch (60.96 cm) seat and a weight bearing capacity up to 400 pounds (204.12

Kg) when evenly distributed. Adjust positioning and unique caregiver only stabilization system for comfort and safety. Tente 5 inch (12.70 cm) casters: three locking and one directional caster. A unique Ultra-Flex suspension system in the seat and back rest cradles the resident from head to heel, greatly reducing the risk of "bottoming out". Ultra-Flex Suspension Bands can be adjusted for individual resident comfort by simply removing the back upholstery. Extra-wide padded armrests and side panels and quick-release upholstery for easy cleaning and replacement. Meets California Technical Bulletin 117-2013. See Swatch Card for available colors/fabrics. Made in the USA.

Item No.: 2959M1



Accessories

**Headrest Pillow:** E-mail [sales@medicus-health.com](mailto:sales@medicus-health.com) for details and pricing

**Adjustable Footrests:** E-mail [sales@medicus-health.com](mailto:sales@medicus-health.com) for details and pricing

**IV Pole:** E-mail [sales@medicus-health.com](mailto:sales@medicus-health.com) for details and pricing

	inch	cm
Seat Width	24.00	60.96
Seat Depth	19.50	49.53
	lb	Kg
Weight Capacity EVENLY DISTRIBUTED	400	181.44

Position  
Stabilization System

Caregiver  
Actuated Position  
Stabilization System

