

Family Practice Table

MODEL 2100M6

- Premium, vinyl upholstery over 2" firm, foam padding (5 cm)
- Full 28.5" width for patient comfort (72.39 cm)
- Air-spring adjustable backrest
- Coated, two position adjustable stirrups
- Pull-out leg rest on steel track
- Easy care all laminate base
- Color coordinated doors and top available
- Chip resistant edges
- 1 large side storage compartment with shelf
- Storage compartment is accessible from both sides (no right or left handed tables)
- 2 drawers with all steel slides
- Built with plywood and laminate panels that are California
- CARB II compliant
- Top sections are easily replaceable in the field
- Adjustable feet for leveling
- Paper dispenser included
- 400 lbs. load capacity under normal use (181.43 kg)



SPECIFICATIONS

58"/72" extended L x 28.5" W x 31" H



California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

Features




Pull-out, slate, laminate, leg rest



2 drawer storage



2-Position Stirrups



Paper Roll Storage



Airlift Backrest



Both Sides Access

Laminate Colors



Dark Cherry



Maple




Gray



Slate Gray

Option



040- Paper Cutter

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black